

Moticam X5 PLUS

WI-FI CAMERAS

Turn your microscope into a Wi-Fi station and share your images with up to 6 devices at the same time. Simply logon with your tablet or smartphone and get ready to capture, edit and store your microscope images with the free MotiConnect App.

The Moticam X5 Plus transfers live images up to 30 fps and captures images at 4.0MP. It also allows you to integrate it to your own network either via Wi-Fi or Ethernet cable, allowing more users to be connected. Alternatively, it also works with our standard Motic Images Plus Software for computers.

<https://www.motic.com/en/>

STREAMING WI-FI CAMERA





A&P INSTRUMENT CO., LTD.
科 藝 儀 器 有 限 公 司

地址：香港九龍灣啟祥道 9 號信和工商中心 1 樓 68 室

Room 68, 1/F Sino Industrial Plaza,

9 Kai Cheung Road, Kowloon Bay, Kowloon

Tel: (852) 2755 6578

Fax: (852) 2755 4549

E-mail: info@anp.com.hk

<http://www.anp.com.hk/>




MOTICONNECT & MOTIC IMAGES PLUS 3.1 SOFTWARE




MotiConnect

View, edit, capture and store your images in the blink of an eye, view live images on multiple wireless devices and share your work with anyone in just a few clicks. The MotiConnect App makes data sharing possible in a collaborative environment, encouraging teamwork and collective learning experiences.



THE [AllinOneBOX] CONCEPT

/ Everything that you may need to work with your Moticam is included in one box



Motic Images Plus 3.1

By connecting your Moticam to your computer you will be able to work with our well-known software, that comes included in the package. View, capture, edit, measure, make reports... all its standard features and the new ones packed in a new user-friendly interface.

Moticam X5 PLUS TECHNICAL SPECIFICATIONS

Sensor	1/3" CMOS Sensor
Imaging Area	5.44 x 3.07mm
Capture Resolution	4MP
Pixel Size	2.0 x 2.0 μ m
Scan mode	Progressive
Shutter mode	Rolling Shutter
Frames per second (fps*)	30fps depending on Wi-Fi strength
Data transfer	Wi-Fi (2.4 or 5 GHz) and Ethernet (RJ45)
Operating temperature	From -10 to +50 Degrees Celsius non condensing
Operating humidity	30~80% RH
Lens mount	CS-Mount
USB	USB 5V Output
RJ45	LAN Port, Default IP: 192.168.2.151
WiFi	2.4G/5G (Support 802.11ac), Default IP: 192.168.1.151
Image Compression Engine	H264 or MJPEG
Stream Resolution	Major: 1920 x 1080, 1280 x 720 Minor: 640 x 360
Software	MotiConnect for Android/iOS/Windows/Mac OS Motic Images Plus 3.1 Windows/Mac OS/Linux
Package includes	12V Power Adapter, CS Ring, Motic 4-dot calibration slide, Dust Cap, Manual, 12mm lens, Adapter φ30/φ38, Macro Tube, Macro Calibration card

*frames per second under optimal illumination conditions